

KaVo PROPHYpearls®.  
Cleaning powder.




**Sodium free!**

New benefits in tooth cleanliness.



KaVo. Dental Excellence.



## A matter of form.

The sphere is the geometrically perfect structure. No other shape has an infinite number of points all the same distance away from the centre.

So the spherical PROPHY-pearls have infinite contact with the contaminated tooth. And they are infinitely gentle in the buccal cavity.

## A clear matter.

With the new PROPHYpearls, tooth cleaning runs true to form. The PROPHYpearls based on lump-free calcium remove deposits, coatings and tartar by their rolling motion - accurately, gently and powerfully. So PROPHYpearls are quickly profitable.

## A rounded matter.

In contrast to conventional sharp-edged cleaning particles, PROPHYpearls are smooth, infinitely soft and gentle to the tooth and its surrounding tissue.



Conventional sodium bicarbonate cleaning powder.




PROPHYpearls without corners and edges.

Order your free PROPHYpearls sample at [www.kavo.com/prophypearls](http://www.kavo.com/prophypearls)


## A clean matter.

Extremely efficient, PROPHYpearls roll over the tooth with an infinite number of contact points - and roll and roll and roll ...



## A matter of fact.

Sharp particles strike the tooth at right angles, bounce off immediately after a single impact. PROPHYpearls, on the other hand, are released at an acute angle for full cleaning power and efficiency.



Sodium bicarbonate:  
90 degree working angle.

PROPHYpearls:  
10 to 60 degree working angle.

DSY DS-No. 3441/11/03 e Subject to alterations.



KaVo. Dental Excellence.

KaVo. D-88400 Biberach/RiB  
Telefon +49 73 51 56-0 · Fax +49 73 51 56-14 88  
Internet: [www.kavo.com](http://www.kavo.com)