

ORIGINAL AIR-FLOW® METHOD
PRODUCT LINE

ORIGINAL AIR-FLOW® METHOD → COMPLETE PRODUCT LINE FOR SUCCESSFUL PROPHYLAXIS

FROM CLASSIC APPLICATION
TO SUBGINGIVAL →
ALL FROM THE NO 1 IN PROPHYLAXIS

ORIGINAL AIR-FLOW® METHOD
SUBGINGIVAL

AIR-FLOW MASTER® THE WORLD'S FIRST SUBGINGIVAL PROPHYLAXIS UNIT

- > For sub- and supragingival treatments – two application systems in one
- > For easy and hygienic operation – with unique 3-touch panel



reddot design award
winner 2008

GOOD DESIGN AWARD 2007
The Chicago Athenaeum
Museum of Architecture and Design



TOUCH
'N'
FLOW

AIR-FLOW® handy PERIO FIRST SUBGINGIVAL PROPHYLAXIS DEVICE FOR THE TREATMENT UNIT

- > For subgingival biofilm management with unique PERIO-FLOW spray handpiece and single-use PERIO-FLOW nozzle



ORIGINAL AIR-FLOW® METHOD
SUPRAGINGIVAL

AIR-FLOW® handy 2+ OUTSTANDING DESIGNER HANDPIECE

- > For more flexibility when using the treatment unit



reddot design award
winner 2003

GOOD DESIGN AWARD 2004
The Chicago Athenaeum
Museum of Architecture and Design



AIR-FLOW® S1 CLASSIC PROPHYLAXIS UNIT

- > For dental plaque removal and tooth polishing



AIR-FLOW® CART MOBILE POWER UNIT

- > For the supply of power and air – with embedded tubing

Excellence
Design Center Stuttgart

Focus mobility 2001



AIR-FLOW® S2 COMBINED PROPHYLAXIS UNIT

- > For removing plaque and for the removal of calculus with the Original Piezon® Method



ORIGINAL AIR-FLOW® METHOD
PROPHYLAXIS POWDERS

AIR-FLOW® POWDER PERIO

SUBGINGIVAL AIR-POLISHING POWDER

- > For the subgingival use of the Original AIR-FLOW® Method
- > For reducing bacteria even in deep periodontal pockets – reduces pocket depth



- > Grain size ~ 25µm (d50)
- > Does not abrade dentin
- > Extrafine

120g bottle

AIR-FLOW® POWDER SOFT

FINE AIR-POLISHING POWDER

- > For the extragente use of the Original AIR-FLOW® Method with regular prophylaxis and recall
- > For use on patients with sensitive gums



- > Glycine-based
- > Grain size ~ 65µm (d50)
- > Does not abrade enamel
- > Extragente

200g bottle

AIR-FLOW® POWDER CLASSIC

ORIGINAL AIR-FLOW® POWDER WITH SIX DIFFERENT FLAVORS

- > For the classic use of the Original AIR-FLOW® Method with basic professional tooth cleaning
- > For removing tooth deposits and stains



- > Sodium bicarbonate
- > Grain size ~ 65µm (d50)
- > Minimally abrasive
- > Anti-inflammatory action

300g bottles

© EMS SA FA-323/EN Edition 04/2012

"I FEEL GOOD"

EMS-SWISSQUALITY.COM

EMS
ELECTRO MEDICAL SYSTEMS SA
Chemin de la Vuarpillière 31
CH-1260 Nyon

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-dent.com

EMS 
ELECTRO MEDICAL SYSTEMS