

# AIR-FLOW®

# 3.0



AIR-FLOW® handy 3.0

**EMS**   
ELECTRO MEDICAL SYSTEMS

AIR-FLOW® handy 3.0

ORIGINAL AIR-FLOW® THERAPY

# AIR-FLOW® handy

ORIGINAL AIR-FLOW® THERAPY WITH A  
AND FLEXIBILITY, ERGONOMICS AND DES



# 3.0

A NEW STANDARD IN EFFICIENCY  
DESIGN



AIR-FLOW<sup>®</sup> PLUS

## PREMIUM PROPHYLAXIS 3.0

MOBILE PROPHYLAXIS DEVICE MEETING THE MOST  
STRINGENT REQUIREMENTS OF MODERN  
PROFESSIONAL TOOTH CLEANING OVER THE ENTIRE  
TREATMENT SPECTRUM BOTH **SUB- AND SUPRAGINGIVAL**

# MODERN PROPHYLAXIS

**AIR-FLOW® handy 3.0**


ORIGINAL AIR-FLOW® THERAPY

## NEW OPTIONS FOR SUCCESSFUL THERAPY →



### **AIR-FLOW® handy 3.0 PREMIUM**

> AIR-FLOW handy 3.0 PLUS with additional AIR-FLOW PERIO handpiece – always ready for sub- and supragingival use with the AIR-FLOW® Powders PERIO and PLUS



### **AIR-FLOW® handy 3.0 PLUS**

> Optimal for use in periodontal pockets up to 4 mm in depth as well as over the entire visible tooth area – flexible and efficient with the AIR-FLOW® Powder PLUS



### **AIR-FLOW® handy 3.0 PERIO**

> For optimal subgingival use in pockets deeper than 4 mm – simple and thorough removal of biofilm with the AIR-FLOW® Powders PERIO and PLUS



### **AIR-FLOW® handy 3.0**

> Specialized in all supragingival applications – for optimal results on the entire visible tooth area with the AIR-FLOW® Powders CLASSIC, CLASSIC COMFORT and SOFT



# BEAUTIFUL → BEAUTIFULLY LIGHT

DESIGN AND FUNCTIONALITY  
FOR OPTIMUM  
TREATMENT COMFORT

- > Designed to ensure easy vision of the treatment area




- > Well-balanced for comfortable handling



## AIR-FLOW® handy 3.0

**AIR-FLOW® handy 3.0**


ORIGINAL AIR-FLOW® THERAPY



**ORIGINAL AIR-FLOW®  
PERIO THERAPY  
PENETRATES INTO THE  
DEEPEST PERIODONTAL  
POCKETS – ELIMINATES  
BIOFILM AND WITH IT  
MILLIONS OF POTENTIAL  
PATHOGENIC GERMS**

PERIO nozzle for subgingival application – three horizontal nozzle outlets for powder-air mixture, one vertical nozzle outlet for water

## **DENTAL IMPLANTS → PROPHYLAXIS WITH THE ORIGINAL AIR-FLOW® THERAPY AND THE AIR-FLOW® handy 3.0 PERIO**




> Implant exposed due to bone loss

- > Peri-implantitis is the periodontitis of implant patients
- > Plaque adheres to the implant and causes inflammation around the implant, which in turn causes bone loss and loosening of the implant
- > Regular subgingival prophylaxis using the Original AIR-FLOW® Therapy prevents peri-implantitis



**AIR-FLOW®  
EASY-FILL**  
> Safe,  
quick and  
clean



## FULL-MOUTH AIR-FLOW® THERAPY

A MIX OF AIR AND ORIGINAL AIR-FLOW®  
POWDER\* COMBINES WITH A JET OF WATER FOR  
RAPID AND PAINLESS REMOVAL OF DENTAL PLAQUE  
AS WELL AS BIOFILM FROM THE ENTIRE MOUTH –  
MINIMALLY ABRASIVE AND DUST FREE

\* All Original AIR-FLOW® Prophylaxis Powders are developed in EMS internal laboratories and are produced by EMS – for optimal performance with Original AIR-FLOW® handpieces



**EMS-SWISSQUALITY.COM**

EMS  
ELECTRO MEDICAL SYSTEMS SA  
Chemin de la Vuarpillière 31  
CH-1260 Nyon

Tel. +41 22 99 44 700  
Fax +41 22 99 44 701  
welcome@ems-ch.com  
www.ems-dent.com