


KaVo. Dental Excellence.

KaVo PROPHYpearls® – perfection all around - only with spheres.

Tried and tested for 10 years.







So comfortable for your patients:

the neutral taste of PROPHYpearls.

Scientifically proven - No adverse modification of tooth t

The Institute for Occupational, Social and Environmental Me PHYpearls cause any scratches or grooves on the blasted dentine or the enamel surface. In processed enamel surfaces, existing grooves



Smooth pearls, smooth result!



KaVo PROPHYpearls -

- Thorough and gentle thanks to the perfect spherical shape
- Safe and bio-compatible due to the calcium carbonate base
- Noticeably pleasant for your patients
- Simple and efficient in their application
- Flexible with the flavours

That brings the fun to prophylaxis!

KaVo's PROPHYpearls all advantages at a glance:

- and has a positive effect on oral flora
- Gentle on the gingiva thanks to
- Pleasant fruity and varied flavours:





Professional prophylaxis with KaVo – unique, versatile and exclusive.

That brings the fun to prophylaxis.

More and more patients realise the importance of reqular prophylaxis, to achieve healthy teeth and a bright smile. Above all, patients expect a pleasant and gentle treatment. KaVo, as a comprehensive service provider for professional prophylaxis solutions, supports you with first-class products for your broad prophylaxis offer.

All of KaVo's prophylaxis solutions inspire through their:

- Simple use
- Well-received treatments for your patients
- Perfect integration into your workflow

Prophylaxis with KaVo simple for you and pleasant for your patients.



10 years KaVo PROPHYpearls now new in many fresh flavours.

For 10 years, KaVo's PROPHYpearls have been proven themselves as reliable and gentle on cleaning powder. The smooth surface of the pearls is particularly convincing. It provides for an effective and gentle tooth cleaning. Now we have improved the KaVo PROPHYpearls even further and are able to offer you the best KaVo PROPHYpearls ever - in addition to many yummy

Whether supragingival or subgingival - KaVo always has a suitable solution.

For thorough and gentle subgingival treatments, KaVo offers its PROPHYflex Perio Powder. This fine-grain glycine powder is particularly gentle on the gingiva and therefore suited best for treatment of periodontitis or peri-implantitis sufferers.

KaVo's range of prophylaxis powders is ideally suited for a multitude of applications:


- Removal of discolorations, plaque and biofilm
- Dental cleaning prior to sealing of fissures
- Dental cleaning prior to bleaching
- For prosthetic and orthodontic applications.
- · For subgingival cleaning and prevention of periodontitis and peri-implantitis





KaVo's range of prophylaxis powders all advantages at a glance:

- An innovative range of powders in new flavours
- For supra and subgingival applications
- Exclusive: the new KaVo PROPHYpearls orange, peach, mint and black currant



KaVo's PROPHYflex powder consists of sodium bicarbonate, with a particle size of between 60μ to 70μ and represents the classic of cleaning agents. Plaque may be easily and efficiently removed while at the same time offering a fruity flavour. KaVo's PROPHYflex powder comes in cherry, berry, mint and orange flavours.













KaVo PROPHYflex® Perio Powder – for gentle thoroughness.

Treat periodontitis and peri-implantitis successfully.

With KaVo's PROPHYflex Perio Powder you can now also work in subgingival regions particularly thorough and gentle. This way you can treat periodontitis and peri-implantitis with long-term success.

The fine-grain glycine powder is particularly friendly to the gingiva and a treatment is significantly more pleasant and better tolerated as e.g. treatment by curette.

The application and handling of the KaVo PROPHYflex Perio Powder is simple and flexible - it may be used with KaVo PROPHYflex 3 2018 or other commercially available powder blasters.



KaVo's PROPHYflex Perio Powder with particularly fine grain is biodegradable (grain size 18–22 μ).

KaVo's PROPHYflex Perio Powder all advantages at a glance:

- horough: clean, efficient cleaning of pockets through removal biofilm (periodontitis-peri-implantitis), reduction of bacterial counts and reduction of the depth of periodontal pockets
- Gentle: finer powder (18-22µm), glycine (active ingredient basis), less abrasive and neutral flavour
- Economic: KaVo PROPHYflex Perio Powder may be used with any commercially available powder blaster. the practical supply bottle is ideal for multiple applications



KaVo's PROPHYflex® powder – where cleanliness has a fruity taste.


Easy and flexible for numerous indications.

With KaVo PROPHYflex powder, you can clean fissures particularly gently prior to sealing; in order to prevent dental caries and gingivitis, you simply and effectively remove plaque, bio-film and discolorations. KaVo PROPHYflex powder is also perfectly suited for use in aesthetic prothetics, for matching the colours of healthy teeth and bridges, inlays or crowns as well as prior to bleaching.









KaVo's PROPHYflex Powder all advantages at a glance:







PROPHYflex powder® suits the taste of your patients.







Fun flavour combinations!



	KaVo PROPHYpearls calcium carbonate	Order number	
NEW!	neutral 4 bottles x 8.8 oz. (250 g)	1.010.1798	
NEW!	neutral 80 sticks x 0.5 oz. (15 g)	1.010.1826	4
NEW!	mint 80 sticks x 0.5 oz. (15 g)	1.010.1828	
NEW!	peach 80 sticks x 0.5 oz. (15 g)	1.010.1829	PEAR
NEW!	orange 80 sticks x 0.5 oz. (15 g)	1.010.1830	
NEW!	black currant 80 sticks x 0.5 oz. (15 g)	1.010.1831	
	KaVo PROPHYflex Pulver Natriumbicarbonat	Order number	0
	KaVo PROPHYflex powder orange 80 Sticks x 0.5 oz. (15 g)	1.007.0014	LVER
	KaVo PROPHYflex powder berry 80 Sticks x 0.5 oz. (15 g)	1.007.0015	
	KaVo PROPHYflex powder cherry 80 Sticks x 0.5 oz. (15 g)	1.007.0016	
	KaVo PROPHYflex powder mint 80 Sticks x 0.5 oz. (15 g)	1.007.0017	
	KaVo PROPHYflex Perio Powder (Glycine)	Order number	
NEW!	KaVo PROPHYflex Perio Powder 4 bottles x 3.5 oz. (100 g)	1.009.3732	

Order number		
1.006.9926	K	
1.006.9927		
1.000.4672	A service (
1.006.9920		
	1.006.9926 1.006.9927 1.000.4672	





