

# PATIENT CHAIR "CECILIA"

HEALTH CARE



MADE IN SWEDEN

## PATIENT CHAIRS "CECILIA"

- Ergonomic quality chairs for efficient examination and treatment
- Easy operation to adjust seat and back support
- New and safe armrest for sampling; "QuickTight"
- Adjustable headrest for ear treatment

The patient chairs called "Cecilia" have adjustable armrests. Height and angle of seat and back support are also adjustable with a simple hand grip.

The Cecilia chairs are available in different designs. When used for examination they are equipped with Rini's new armrest for sampling; QuickTight. They are easy to adjust and lock in desired position, comfortable for both patient and personnel.

Cecilia is an ear examination chair equipped with a headrest which is adjustable in both height and depth.

The chairs are provided with trumpet gliders and dressed in beige PVC fabric. As an alternative to trumpet glider, the chairs can be supplied with wheels.



## TECHNICAL INFORMATION

### GENERAL

Length (transport mode)	580 mm
Width (transport mode)	530 mm
Weight unit	20 kg/44 lbs
Height adjustment from floor	450-590 mm
Lifting capacity	120 kg/265 lbs

### BACK SUPPORT

	"Cecilia Test"	"Cecilia Ear"
Height	230 mm	230/130 mm
Width	380 mm	380/260 mm
Height adjustment from seat	40-150 mm	40-150 mm

### SEAT

Length	420 mm
Width	460 mm



Rini Ergoteknik AB Fågelsångsvägen 3B  
SE-186 42 Vallentuna Sweden  
Tel +46 8-594 771 70 Fax +46 8-929 279  
info@rini.se www.rini.se