

# Medical-D

Medical-D digital ring flash camera  
for easy dental imaging

trademark

## 3 dial modes


SMILE mode  
D-mode (close-up)  
F-mode (close-up +  
anti-glare)  
> use optional CP filter)

## ring flash light



## MACRO ring flash

detachable  
distance positioner



SONY DSC-HX99  
18.2 Megapixels  
Close-up Macro +2  
CMOS sensor

- quick and Easy
- WiFi certified
- Micro SD XC
- automatic light
- battery charged
- Use add-on CP filter;  
[Separately available]
- practical bag
- touch focus
- tilting rear screen



Medical-D is  
produced by

DHM-dental B.V.  
Simon Stevinweg 48  
6827 BT Arnhem  
the Netherlands

93/42 EEG  
MDR:2017/745



additional polarisation filter available

www.dhm-dental.com

- \* easy to use
- \* Medical-D wifi
- \* full HD film
- \* tilting touch screen
- \* protection bag

#### Included in set



Li-ion



charger



manual



warranty



store-bag

**F3.5**

**18.2MP**  
CMOS

**3.0"**  
LCD

**2x macro**  
close-up

**30X**  
ZOOM



practical  
protection bag



# Medical-D

Anyone can take perfect photos for dental purposes with the Medical-D camera, quickly and simply, by using the special pre-programmed modes (D-mode, smile-mode + optional F-mode), for mouth and / or face images.

If you want to make anti-glare images with the prepared F-mode function you can use the optional 'add-on' crosspolarisation filter (separately available).

It is also possible to record a full HD video of the relationship or articulation. Direct access to history and any picture enlargements by the touch screen on the backside.

This camera uses built-in WiFi which enables photos to be sent to a PC, laptop, tablet or mobile phone directly.

A unique camera because of its speed, simplicity and quality.



**SONY DSC-HX99**  
18.2 Megapixels  
Close-up Macro +2  
30 x ZOOM

Medical-D® is a registered trademark.

This product is designed, developed and produced in the Netherlands by DHM-dental BV.