# Technical data



SURGmatic S15L Pro			
1.014.4000	1.013.7541	1.013.7540	1.009.1010
1:5	20:1	20:1	1:1
40,000 rpm	40,000 rpm	40,000 rpm	40,000 rpm
1,500-200,000 rpm	15-2,000 rpm	15-2,000 rpm	300-40,000 rpm
No	Yes	No	No
N/A	80 Ncm	55 Ncm	5.5 Ncm
Yes	No	No	N/A
Yes	Yes	Yes	Yes
Yes	Yes	Yes	Yes
Yes	Yes	Yes	Yes
No	Yes	Yes	No
Yes	Yes	Yes	Yes

<sup>\*</sup> with KaVo surgical motor S600LED; motor speed range 300-40,000 rpm

# Associated products

The system concept is very important to us: in order to make everyday practice as comfortable for you as possible, we ensure that our products are seamlessly compatible.



#### KaVo MASTERsurg

Convenient operation through the intuitive touchscreen display, the freely positionable wireless foot control and the data documentation feature with a graphic real-time display showing torque for optimal quality control.



#### KaVo EXPERTsurg

Rely on high performance with 80 Ncm torque on the instrument and save time through intuitive programming of your work steps and parameters.

The products, equipment and services shown and described in this catalogue are not available in all countries. All information corresponds to the knowledge available at the time of printing. KaVo Dental GmbH assumes no responsibility for discrepancies in the colour and shape of the images, or for errors or misprints, and reserves the right to make changes to the brochures at any time. Reprinting, even in part, is only authorised with the approval of KaVo Dental GmbH.





# SURGmatic S15 L Pro: The new all-rounder amongst surgical instruments.

## NEW

# At top speed for more efficiency.

## Speed translated into cutting performance

With a ratio of 1:5 and up to 200,000 revolutions per minute, the SURGmatic S15 L Pro sets a new standard in surgical dentistry and takes surgical work to a whole new level.


#### A real all-rounder

Its advanced speed asset makes the KaVo SURGmatic S15 L Pro not only noticeably more efficient, but also a real all-rounder in the surgical practice.

In addition to surgical applications such as the removal of wisdom teeth, hemisections, osteotomies or apicectomy, non-surgical activities such as preparations, crown lengthening or the separation of crowns can also be carried out effortlessly.

Thanks to the high speed and the compelling cutting performance, inconvenient instrument changing is no longer necessary and a smooth workflow is guaranteed.






 $<sup>^{\</sup>star}$  Representation of the average period of time for cutting through a defined sample of zirconium oxide

# Working simply becomes easier - you'll have a cool head and a clear view.

## Hotheads don't stand a chance

The new KaVo CoolHead technology actively prevents excessive heating of the head if the push button is accidentally pressed.





## Compelling ergonomics, hygiene and quality

With the compact head size, the low weight and the compelling lighting effects, the SURGmatic S15 L Pro creates the best conditions for ergonomically healthy and permanently relaxed work.

The precise KaVo Made in Germany manufacturing quality, with high-quality materials and a hard metal guide bushing ensure more precision thanks to low vibrations and a high level of robustness. Of course, the SURGmatic S15 L Pro can be thermodisinfected and sterilised up to 135 °C.



enables short, gentle

patient treatment.

Everything under one click:
See the new SURGmatic S15 L Pro
for yourself online.

www.kavo.com/en/s15