

This table steam generator, compact and practical, is ideal for sanitation. It respects the tradition of sturdiness, practicalness and long-lasting life of our products. Connection to the water mains is not required. Quick and easy to use: just pour decalcified water into the tank and plug the steam generator in; it will be ready in approximately 8 minutes, then all you will have to do is to top up the water tank whenever indicated by the pilot light and the cycle will continue without interruptions. Hands-free equipment, thanks to the foot control, for a practical and quick use.



**Colours available on stock**



*Oro*

*Argento*

*Bianco*

**OPTIONAL: colours available on request**



*Nero*

*Bordeaux*

*Arancio*

*Verde lime*