


**ONE  
FOR  
ALL**



**GUIDED  
BIOFILM  
THERAPY**



## BIOFILM MANAGEMENT WITH AIR-FLOW® PLUS POWDER



ABOVE  
CEMENTO-  
ENAMEL  
JUNCTION

### 1. AIR-FLOW® SUPRAGINGIVAL



CEMENTO-  
ENAMEL  
JUNCTION  
< 4mm

### 2. AIR-FLOW® SUB- AND SUPRAGINGIVAL



CEMENTO-  
ENAMEL  
JUNCTION  
> 4mm

### 3. AIR-FLOW® SUBGINGIVAL

**“I FEEL  
GOOD”**

**EMS-DENTAL.COM**

EMS  
Electro Medical Systems SA  
Chemin de la Vuarpillière 31  
CH-1260 Nyon  
Phone +41 22 99 44 700  
Fax +41 22 99 44 701  
contact@ems-ch.com

W&H Nordic AB  
Tillverkarvägen 6  
S-187 11 Täby  
Sweden

© EMS FA-569/EN - rev. A-01

ALL-IN-ONE  
**N°1**  
PROPHYLAXIS

ORIGINAL  
METHOD

**AIR-FLOW® POWDERS**

**GUIDED  
BIOFILM  
THERAPY**

**EMS+**

## AIR-FLOW® POWDERS

## ORIGINAL THERAPY AIR-FLOW®

### → SUB- AND SUPRAGINGIVAL



#### AIR-FLOW® POWDER PLUS


- SUB- AND SUPRAGINGIVAL – IN ONE STEP
- NO IRRITATION OF THE GUMS – PATIENTS FEEL GOOD AND FREE OF PAIN
- NO MEASURABLE IMPAIRMENT OF THE TOOTH SURFACE, EVEN IN RECALL PROPHYLAXIS
- Erythritol-based
- Average grain size ~14 µm
- 120 g bottle

### CLINICAL APPLICATIONS

- REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS WITH **OPTIMIZED PATIENT COMFORT**
- ELIMINATES SUPERFLUOUS POLISHING WITH POLISHING PASTE – THAT SAVES TIME

- REMOVAL OF BIOFILM FROM POCKETS **DOWN TO 4 MM** ON NATURAL TEETH AND IMPLANTS
- REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS, ON NATURAL TEETH AND IMPLANTS **FROM 4 MM AND DEEPER**
- WITH AIR-FLOW® PERIO HANDPIECE AND SOFT DISPOSABLE NOZZLE WITH DEPTH MARKS

### AIR-FLOW® HANDPIECES



### → SUPRAGINGIVAL



#### AIR-FLOW® POWDER CLASSIC NEW FORMULA

- Sodium bicarbonate
- Average grain size ~40 µm
- 300 g bottle
- Four flavors: lemon, mint, cherry and neutral

- REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS WITH OPTIMIZED **PATIENT COMFORT** ABOVE CEMENTO-ENAMEL JUNCTION

