

KaVo PHYSIO Evo / Evo F

Sitting comfortably:
Ergonomics meet design.



KaVo. Dental Excellence.

Comfort made beautiful.

It's all about your back.

For dentists, it is particularly important to assume a steady, ergonomically healthy and comfortable body posture in any treatment situation. And nothing influences your body posture as much as the chair that you sit on every day.

KaVo is renowned for creating products that optimise both ergonomic design and comfort. In keeping with this, the KaVo PHYSIO Evo and Evo F is precisely designed according to the needs and workflows of dentists and dental assistants.



Ergonomics and comfort in perfect harmony: with the KaVo PHYSIO Evo, you always sit with a healthy, natural posture.

Correct posture is key.

The intelligent chair concept of the KaVo PHYSIO Evo fulfills the highest demands on an ergonomic sitting position. Thanks to independent adjustment functions, the seating area and the backrest may be set individually to your anatomical requirements and personal treatment style. The greatest possible sitting comfort allows for the smoothest workflow.

Ergonomics are a matter of detail.

The backrest and seat cushion of the KaVo PHYSIO Evo are ergonomically shaped. The central ridge in the seat cushion provides stable sitting and prevents any sliding even in case of short, quick rolling motions. The rounded thigh support allows for freedom of motion.

Thanks to the spring force-supported backrest, the back muscles – in particular in the lumbar region – are permanently actively supported and thus continuously protected.

An excellent concept.

The KaVo PHYSIO Evo has been awarded the seal of approval of the German association "Campaign for Healthier Backs" (AGR). This is awarded by an independent, multidisciplinary expert commission according to strict test criteria regarding ergonomics and health.



Innovative design and KaVo quality go hand in hand.

Form and colour sit well together.

The KaVo PHYSIO Evo has a lean, modern construction and offers a large variety of colours. The compact construction and small diameter of the base offer plenty of freedom of movement even within a tight space. The colour range of the KaVo PHYSIO Evo matches the colours of the treatment units KaVo ESTETICA E80, E70 and E50.

KaVo quality you can count on.

High-quality materials produced to the highest KaVo standards plus comprehensive endurance tests guarantee a long service life. Likewise, the KaVo PHYSIO Evo allows for quick and reliable hygiene: all the components used in manufacture have been assembled with the lowest gaps and the surface is easily cleaned.

The KaVo PHYSIO Evo F is also exceptionally hygienic. Cross contamination is prevented due to the seat height adjustment via foot control.

Dental excellence pays off.

An ergonomically healthy sitting posture may have a significant influence on the profitability of your practice. Because in the long term, continuously working with a healthy posture and not experiencing any health issues is priceless for you and your employees.

All the benefits at a glance.

- Perfect ergonomics with maximum comfort for healthy and fatigue-free working
- Perfectly finished KaVo quality for trouble-free performance in the day-to-day operations of your practice
- Attractive colours and modern design recognised by the GOOD DESIGN AWARD
- Best hygiene thanks to precisely tuned components with very small gaps and surfaces that are easy to clean.



An extremely lean base for plenty of freedom of motion and soft, braked castors for controlled rolling motions.



Technical Specifications

	PHYSIO Evo	PHYSIO Evo F
Minimum height of the seat	465 mm	475 mm
Maximum height of the seat	660 mm	660 mm
Adjustment range of the seat height	195 mm	185 mm
Width of the seat	max. 400 mm	max. 400 mm
Depth of the seat	380 mm	380 mm
Adjustment range for the backrest/seat tilt	30°	30°
Seat inclination	15°	15°
Base diameter	535 mm	535 mm
Weight	12,5 kg	13 kg



We reserve the right to make technical modifications. Slight colour differences are due to the printing process. © Copyright KaVo Dental GmbH.



KaVo. Dental Excellence.