

**NEW!**

**FUN,  
FRESH AND  
EFFECTIVE.**



KaVo. Dental Excellence.

# KaVo PROPHYpearls® – perfection all around – only with spheres.



## Tried and tested for 10 years.

KaVo PROPHYpearls are spherical and thus offer optimal cleaning performance: the round surfaces possess more contact points than conventional powder and thus work even more efficiently. At the same time they are particularly gentle on the teeth, as they do not have any sharp edges. Treatments are very well received by patients. KaVo's PROPHYpearls come in the following fresh flavours: black currant, mint, orange, peach and neutral.

There are many applications for KaVo's PROPHYpearls:

- Removal of discolorations, plaque and biofilm
- Dental cleaning prior to sealing of fissures
- Dental cleaning prior to bleaching
- Prosthetic and orthodontic applications.



KaVo's PROPHYpearls without corners or edges. (Grain size 60–70 µm).



Cleaning by means of KaVo's PROPHYpearls is effective and friendly to the substance. (x150 magnification).



With their infinite number of contact points, KaVo's PROPHYpearls roll across tooth surfaces and therefore work even more efficiently.



KaVo's PROPHYpearls and KaVo's PROPHY Superpearls: angle of operation 10° to 60°. Thanks to this more acute angle, KaVo's PROPHYpearls roll with an infinite number of contact points across the tooth surface and thus release their full effect.

## Cleaning more gently with PROPHYpearls®.

### Perfect form and function:

PROPHYpearls are based on non-clumping calcium and may be used at any time. They are ball-shaped and provide optimum cleaning power: thanks to an infinite number of contact-points on the spherical surface. The dynamic PROPHYpearls will roll across tooth surfaces and will correspondingly work more efficiently.

### So comfortable for your patients:

PROPHYpearls are particularly gentle on the teeth, as they do not have any sharp edges. And your patients will perceive it as extremely pleasant that the tooth cleaning with PROPHYpearls is gentle to the dental environment as well. Furthermore, many patients prefer the neutral taste of PROPHYpearls.

### Save valuable time with PROPHYpearls:

Apart from the actual tooth cleaning, handling the PROPHYpearls is very efficient for you. Simply refill the pearls with the pre-dosed sachet or with the new bottles.

### Scientifically proven – No adverse modification of tooth

The Institute for Occupational, Social and Environmental Medicine in Ulm (Germany) confirms: "There were no indications that PROPHYpearls cause any scratches or grooves on the blasted dentine or the enamel surface. In processed enamel surfaces, existing grooves were not levelled out in any noticeable way." You can find further information at: [www.kavo.com/studien](http://www.kavo.com/studien)

# PEARLS



# Smooth pearls, smooth result!



KaVo PROPHYpearls – the original from KaVo for a perfect and pleasant clean.

Now with new flavours.

- **Thorough** and gentle thanks to the perfect spherical shape
- **Safe** and bio-compatible due to the calcium carbonate base
- **Noticeably pleasant** for your patients
- **Simple and efficient** in their application
- **Flexible** with the flavours  
mint, peach, orange, blackcurrant

That brings the fun to prophylaxis!

KaVo's PROPHYpearls – all advantages at a glance:

- Increases the pH, reduces acids and has a positive effect on oral flora
- Gentle on the gingiva thanks to Microsphere technology
- Pleasant fruity and varied flavours: orange, mint, peach and black currant



**NEW!**





# Professional prophylaxis with KaVo – unique, versatile and exclusive.

The original  
by KaVo.  
**10 Years**  
PROPHYpearls

collection

## That brings the fun to prophylaxis.

More and more patients realise the importance of regular prophylaxis, to achieve healthy teeth and a bright smile. Above all, patients expect a pleasant and gentle treatment. KaVo, as a comprehensive service provider for professional prophylaxis solutions, supports you with first-class products for your broad prophylaxis offer.

All of KaVo's prophylaxis solutions inspire through their:

- Simple use
- Well-received treatments for your patients
- Perfect integration into your workflow

## Prophylaxis with KaVo – simple for you and pleasant for your patients.

## 10 years KaVo PROPHYpearls – now new in many fresh flavours.

For 10 years, KaVo's PROPHYpearls have been proven themselves as reliable and gentle on cleaning powder. The smooth surface of the pearls is particularly convincing. It provides for an effective and gentle tooth cleaning. Now we have improved the KaVo PROPHYpearls even further and are able to offer you the best KaVo PROPHYpearls ever – in addition to many yummy flavours.

## Whether supragingival or subgingival – KaVo always has a suitable solution.

For thorough and gentle subgingival treatments, KaVo offers its PROPHYflex Perio Powder. This fine-grain glycine powder is particularly gentle on the gingiva and therefore suited best for treatment of periodontitis or peri-implantitis sufferers.

KaVo's range of prophylaxis powders is ideally suited for a multitude of applications:

- Removal of discolorations, plaque and biofilm
- Dental cleaning prior to sealing of fissures
- Dental cleaning prior to bleaching
- For prosthetic and orthodontic applications.
- For subgingival cleaning and prevention of periodontitis and peri-implantitis



KaVo's range of prophylaxis powders – all advantages at a glance:

- An innovative range of powders in new flavours
- For supra and subgingival applications
- Exclusive: the new KaVo PROPHYpearls orange, peach, mint and black currant

## The classic with sodium bicarbonate – clean in any case.

KaVo's PROPHYflex powder consists of sodium bicarbonate, with a particle size of between 60µm to 70µm and represents the classic of cleaning agents. Plaque may be easily and efficiently removed while at the same time offering a fruity flavour. KaVo's PROPHYflex powder comes in cherry, berry, mint and orange flavours.





# KaVo PROPHYflex® Perio Powder – for gentle thoroughness.

## Treat periodontitis and peri-implantitis successfully.

With KaVo's PROPHYflex Perio Powder you can now also work in subgingival regions particularly thorough and gentle. This way you can treat periodontitis and peri-implantitis with long-term success.

The fine-grain glycine powder is particularly friendly to the gingiva and a treatment is significantly more pleasant and better tolerated as e.g. treatment by curette.

The application and handling of the KaVo PROPHYflex Perio Powder is simple and flexible – it may be used with KaVo PROPHYflex 3 2018 or other commercially available powder blasters.



KaVo's PROPHYflex Perio Powder with particularly fine grain is biodegradable (grain size 18–22 µ).

### KaVo's PROPHYflex Perio Powder – all advantages at a glance:

- **thorough:** clean, efficient cleaning of pockets through removal biofilm (periodontitis-peri-implantitis), reduction of bacterial counts and reduction of the depth of periodontal pockets
- **Gentle:** finer powder (18–22µm), glycine (active ingredient basis), less abrasive and neutral flavour
- **Economic:** KaVo PROPHYflex Perio Powder may be used with any commercially available powder blaster, the practical supply bottle is ideal for multiple applications

PERIO



# KaVo's PROPHYflex® powder – where cleanliness has a fruity taste.

## Easy and flexible for numerous indications.

With KaVo PROPHYflex powder, you can clean fissures particularly gently prior to sealing; in order to prevent dental caries and gingivitis, you simply and effectively remove plaque, bio-film and discolorations. KaVo PROPHYflex powder is also perfectly suited for use in aesthetic prosthetics, for matching the colours of healthy teeth and bridges, inlays or crowns as well as prior to bleaching.



Prior to cleaning.



Post cleaning with PROPHYflex.



Scanning Electron Microscope image post cleaning with KaVo's PROPHYflex 3 (x20 and x200 magnification).



Conventional sodium bicarbonate: angle of operation 60° to 90°. The sharp-edged particles of the powder impact with the tooth surface and immediately bounce off.

## PROPHYflex powder® suits the taste of your patients.

Respond to the individual taste of your patients. This is possible due to PROPHYflex powder being available in the fruity flavours cherry, berry, mint and orange. KaVo's PROPHYflex powder consists of sodium bicarbonate, with a particle size of 60µ – 70µ. Plaque may be easily and efficiently removed while at the same time offering a fruity flavour.

### KaVo's PROPHYflex Powder – all advantages at a glance:

- Thorough and gentle thanks to fine-grain powder
- Noticeably pleasant for your patients
- Simple and efficient use
- Flexible through flavours orange, berry, mint and cherry.



PERIO



# Fun flavour combinations!

**NEW!**

**NEW!**

**NEW!**

**NEW!**

**NEW!**

**NEW!**

KaVo PROPHYpearls calcium carbonate	Order number
neutral 4 bottles x 8.8 oz. (250 g)	1.010.1798
neutral 80 sticks x 0.5 oz. (15 g)	1.010.1826
mint 80 sticks x 0.5 oz. (15 g)	1.010.1828
peach 80 sticks x 0.5 oz. (15 g)	1.010.1829
orange 80 sticks x 0.5 oz. (15 g)	1.010.1830
black currant 80 sticks x 0.5 oz. (15 g)	1.010.1831
KaVo PROPHYflex Pulver Natriumbicarbonat	Order number
KaVo PROPHYflex powder orange 80 Sticks x 0.5 oz. (15 g)	1.007.0014
KaVo PROPHYflex powder berry 80 Sticks x 0.5 oz. (15 g)	1.007.0015
KaVo PROPHYflex powder cherry 80 Sticks x 0.5 oz. (15 g)	1.007.0016
KaVo PROPHYflex powder mint 80 Sticks x 0.5 oz. (15 g)	1.007.0017
KaVo PROPHYflex Perio Powder (Glycine)	Order number
KaVo PROPHYflex Perio Powder 4 bottles x 3.5 oz. (100 g)	1.009.3732



**NEW!**

PROPHYflex colour	Order number
KaVo PROPHYflex 3 2018 violet	1.006.9926
KaVo PROPHYflex 3 2018 violet Sirona	1.006.9927
KaVo PROPHYflex 3 2018 black	1.000.4672
KaVo PROPHYflex 3 2018 black Sirona	1.006.9920



**KaVo. Dental Excellence.**