New seal²





A perfect seal in 2 seconds!

Seal² is an automatic pouch sealing device that produces 12 mm seals. The sealing temperature is electronically regulated to guarantee a fully hermetic seal and to prevent bags from being burnt.

Comfort and ease of use

Handling is improved with the new transparent sealing station and integrated work surface.

All types of roll in arm's reach

Seal^{2's} patented double roll holder can accommodate a large selection of roll widths to suit all your needs.

Seal² fits everywhere

The roll holder can be separated from the sealing station and can also be wall mounted.



W&H Sterilization S.r.I. Italy, I-24060 Brusaporto Via Bolgara 2 wh.com