

# W&H Couplings

Versatile system interfaces



**Flexible** – for use with all standardised turbine tubing. **Full compatibility** – compatible with all turbines, air motors and scalers with Roto Quick and Multiflex®\* connections. **360° instrument rotation** – for easy handling, prevents the tube from twisting. **Continuous regulation of the spray water supply** – quick, customised adjustment by the user even –during treatment. **Non-retraction valve** – avoids back-flow of contaminated spray water to minimise the risk of infection.



Scan the QR code with your mobile phone for more information on Couplings

Diagrams, symbol photos, additional equipment and the content of illustrated accessories are not included in delivery.



Roto Quick couplings:	RQ-24	RQ-34	RQ-14	RQ-03	RQ-04
Light:	yes	yes	–	–	–
Spray control:	–	yes	yes	–	–
Connection:	6-hole	6-hole	Standard 4-hole	2 (3)-hole	Standard 4-hole
					

### W&H Roto Quick couplings

Light, flexible and sterilizable: W&H Roto Quick couplings are versatile to use and are suitable both for instruments with LED+ and LED, as well as halogen light. 24 months of warranty.

W&H coupling for Multiflex®* connection:	RM-34 LED
Light:	LED
Spray control:	yes
Connection:	6-hole
	



### W&H coupling RM-34 LED for Multiflex®\* connection

Simple and uncompromising: The W&H coupling RM-34 LED is ready for immediate use. This means: The coupling with LED can be used alongside the conventional halogen couplings without requiring any changes at all in the dental unit. Long service life – 24 months of warranty.

\* Multiflex® is a registered trademark of Kaltenbach & Voigt GmbH, Germany

W&H Dentalwerk Bürmoos GmbH  
 Ignaz-Glaser-Straße 53, Postfach 1  
 5111 Bürmoos, Austria  
 t + 43 6274 6236-0  
 f + 43 6274 6236-55  
 office@wh.com  
 wh.com